	FR	DE	NL	EN
①	ON/OFF	EIN/AUS	AAN/UIT	ON/OFF
②	Sens de rotation	Drehrichtung	Rotatierichting	Direction of rotation
③	3 intensités de massage	3 Massageintensitäten	3 massage intensiteiten	3 massage intensity levels
④	Bouton de chauffe	Heiztaste	Verwarmingsknop	Heat button

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FR

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Warning: Please read this instruction manual carefully before using the instrument and keep it for future use.

A. GENERAL SAFETY INSTRUCTIONS

1. The device is intended exclusively for private and not for commercial use. Only use this device as describe in this instruction manual.
2. Before to use, make sure that the supply voltage corresponds with the voltage marked on the rating label.
3. Do not leave the device unattended when plugged in.
4. Do not handle the unit or plug with wet hands.
5. Unplug from the power outlet after use and allow the appliance to cool down completely before to clean it or to move it.
6. When disconnecting the cable, always pull at the plug; do not pull on the cable itself.

7. When using the device, always pull the power cord completely from its housing to prevent overheating.
8. If the cord is damaged, it must be replaced by the manufacturer or its service agent or a similarly qualified person in order to avoid an electric shock.
9. Never try to repair it by yourself. It must be repaired by the manufacturer or its service agent or a similarly qualified person in order to avoid an electric shock.
10. The power cable must be checked regularly to avoid damage. If the power cable is damaged it must be replaced by a qualified service center to avoid any danger.
11. Always place your device on a dry surface.
12. Do not use outdoors.
13. Do not let the power cord on the edge of a table or work surface, and prevent it comes into contact with hot surfaces.
14. Keep it away from heat sources such as radiators, to avoid deformation of plastic parts.
15. Do not touch the hot surface. Keep the appliance out of reach of children.
16. This unit is not designed to be operated by means of an external timer or a separate remote control system.

B. SPECIAL SAFETY INSTRUCTIONS

1. This device must not be used by children under 3 years of age as they are unable to react in the event of it overheating.
2. This device must not be used by small children (3-8 years) unless a parent or guardian has pre-set the necessary controls, or the child has been given instructions on how to use the heating pad safely.
3. This device can be used by children aged 8 or over, and by persons with physical, sensory or mental disability or with limited experience or knowledge, provided they are under supervision or have received instructions on how to use it.
4. Children must not play with it.
5. Children must not be allowed to clean or maintain this device without supervision.
6. Make sure that you do not fall asleep while using this device!
7. If you are very tired while using it, you can select level 1 for your safety.
8. You must examine the device frequently for signs of wear and tear or damage. If you detect signs of wear or tear or damage, or if the device has been used inappropriately stop to use it.
9. Do not attach safety pins, needles or pins of any kind to this device.
10. This product is not intended for use in hospitals, clinics etc
11. This is not a medical product.
12. Do not use with a helpless person, child or any person with low sensitivity to heat.

C. WARNING

1. Persons with the following conditions should not use this product or should consult their physician before using the product:
 - Pregnant women
 - For those with an integrated medical device, such as a pacemaker.
 - Those with diabetes
 - Suffering from a blood circulation problem or thrombosis.
 - Suffering from a skin disease or infection.
 - Having wounds, burns, sores or other skin alterations.
 2. Do not use on swollen or inflamed areas.
 3. Stop using immediately if you feel discomfort or pain.
 4. Do not use this device before going to bed. Massagers are stimulating and may cause delays in sleep.
 5. To minimize the risk of electric shock, do not use on wet skin.
 6. Keep the device away from moisture. Do not use the device or any of its components if they are wet.
 7. Always keep the equipment clean. Do not apply lotions, creams, etc. to the skin before or during use.

D. USE

- Plug the adapter into a standard electrical outlet.
- Place the massager on the floor, in front of a chair or seat.
 - Press the ① button to turn on the device.
 - Press the ② button to change the direction of rotation.
 - To adjust the intensity, choose one of the 3 modes by pressing the ③ button.
 - Press the heating ④ button to begin heating. The heating indicator light will turn red.
 - Press the ④ button again to turn off the heating function.

M 3200 / Z5

- The device automatically turns off after 15 minutes of use. To continue the massage, press the ① button, the device restarts.

Caution: Be sure to unplug the unit from the power outlet immediately after use and before cleaning. It is strongly recommended that people wear socks when using the foot massager. Do not use the unit with shoes on. The unit has a built-in timer with a safety feature that, after turning off the unit, stops the heating and massage after 15 minutes. It is recommended to wait 15 minutes before turning it back on.

You can also unzip the removable part ⑥ to use the cushion alone ⑦.

E. CARE AND MAINTENANCE

Wipe only with a soft, slightly damp sponge. Never allow water or any other liquid to come into contact with the unit.

Do not use abrasive cleaners to clean the unit.

Store in a cool, dry place out of direct sunlight. Avoid humidity.

F. TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
The device will not operate at all	might not be plug or start on	Plug the adapter into the socket again. Check the operation of the button.
The device shuts down unexpectedly	The device may have been used for 15 minutes. The automatic shutdown function is activated.	Press the "ON" button after turning off the device for 15 minutes to continue using it.

G. FEATURES

- Material: faux leather, mesh
- Power supply : DC 12V (power adapter included)
- Massage: 3 intensities
- Heat: ON/OFF 1 heat intensity level

The European directive 2012/19/EC about electronic and electric waste, requires that you can't throw away defective domestic appliance with common waste. Used device has to be collect separately to optimize the recuperation rate and the recycling of the materials to reduce the impact on the health and the environment.

Importé par / Importiert von / Geïmporteerd door / Imported by :
CAPTELEC, 59170 Croix - FRANCE